

The DISTEL DOSING SYSTEM is intended for the cleaning and disinfection of surfaces of non-invasive medical devices and equipment, and of hard non-porous environmental surfaces.

A 500ml bottle of DISTEL DOSING SYSTEM concentrate generates measured doses of surface disinfectant solution for application via a Doser or Sprayer.

The DISTEL DOSING SYSTEM working solution is effective against coronaviruses in 2 minutes. Efficacy is also proven against bacteria, yeasts, enveloped viruses and nucleic acids (DNA/RNA)*


*Please see the microbiological efficacy summary for specific contact times.

ELIMINATE SINGLE-USE PLASTIC WIPES

Designed to be repurposed and reused, the DISTEL DOSING SYSTEM keeps the environment in mind.

A 500ml bottle of DISTEL DOSING SYSTEM concentrate delivers 10 x 500ml bottles of working solution.

The DISTEL DOSING SYSTEM can be used with any spreader to create freshly made surface disinfectant wipes on the go. Choose a sustainably sourced spreader to reduce single-use plastic waste.



500ML BOTTLE OF CONCENTRATE







ONE SOLUTION TWO APPLICATORS

The DISTEL DOSING SYSTEM puts you in control.

Choose from two applicators:

- The Doser, for measured doses of liquid
- The Sprayer, for measured doses of liquid spray



THE DISTEL DOSING SYSTEM IS INTENDED FOR THE CLEANING AND DISINFECTION OF SURFACES OF NON-INVASIVE MEDICAL DEVICES AND EQUIPMENT, AND OF HARD **NON-POROUS ENVIRONMENTAL SURFACES SUCH AS:**



Training and certification is an essential part of customer service with Cache. Ask about our Online Training Portal for product training at a time and place convenient to you.

AVAILABLE TO ORDER NOW!



6 x 500ml Dosing Bottles

3 x SPRAYER

CCH030501 3 x DOSER

Bottles - 3 x 5ml Doser 3 x 1ml Sprayer - 1 x Working Solution Bottle Label

cache

1 x 500ml Bottle of Concentrate - 5 x 500ml Dosing

SYSTEM CCH050501 1 x 500ml Bottle of Concentrate 1 x Working Solution Bottle Label

